

## Helper cushion



The dental professional posture is dictated by what they need to see.

In order to reduce high risk back and neck flexions, patient head must be at the end of headrest.

**Dimensions** (width, length, height):  
18 x 18 x 6 in (46 x 46 x 15 cm)

**Material:** foam and antibacterial upholstery

- **Better visibility:** Keeps patient close to you
- **Better posture:** Improves professional posture and patient comfort
- **Dental chair complement:** To be used with most dental patient chairs
- Straps for easy storage
- Resistant to disinfection

**Maintenance:** Our regular upholstery (Naugahyde, Naugasoft) contains agents to protect it against bacterial and fungal micro-organisms. A soft soap and water are to be used for regular cleaning. When disinfection is required, use a stabilized hydrogen peroxide disinfectant or a low alcohol concentration product (less than 10 %) for all components.

## Suggested use

